

Usage of a Power Inverter

Sample Power Inverter:



Black and Decker
400watts
Powers to 3.48 AMPS

Note: You may use other brands that are available from local auto parts stores, retail stores and online however, you must purchase at least a 400watt inverter and one that powers to 3.48AMPS or more.

When using a Power Inverter for traveling with the CryoPen in between demo appointments, please make sure to do the following:

1. Attach inverter to vehicle following manufacturer's instructions.
2. Plug CryoPen into inverter
3. Turn car on
4. Turn inverter on
5. Turn CryoPen on
(Make sure the position of the CryoPen when traveling, does not block the air flow to the unit. If airflow is blocked, unit will overheat and become inoperable).

Arriving at destination:

6. Turn CryoPen off and unplug from inverter
7. Turn inverter off
8. Turn car off

Do Not turn car off before turning unit off, unplugging from inverter and turning inverter off.

Turning car off first will render the unit inoperable. Car must be running while turning CryoPen and inverter off.